# SANAO

## Handpiece

### Instructions for use

EN  
DE  
FR

CS / IT / NL / ...

via: [www.scican.uk.com/downloads](http://www.scican.uk.com/downloads)

200 L  
40 / 40 L  
10 / 10 L  
PSO  
PSI  
40 ST

Rx only

# SciCan

Your Infection Control Specialist™

[www.scican.com](http://www.scican.com)

Your Endo Specialist™

[www.micro-mega.com](http://www.micro-mega.com)

# SANAO

Copyright 2016 SycoTec. All rights reserved.

Doc. No. 2.002.1189 - Rev. 2016-04

SANAO and STATMATIC are trademarks of SciCan Ltd.  
All other trademarks referred to in this handbook are the  
property of their respective owner.

Produced by:

---

**Sycotec GmbH & Co. KG**

Wangener Strasse 78  
88299 Leutkirch  
GERMANY  
Telephone: +49-7561-86-0  
Fax: +49-7561-86-400  
E-mail: info@sycotec.eu

Exclusively distributed by:

---

**SciCan GmbH**

Wangener Strasse 78  
88299 Leutkirch, GERMANY  
Telephone: +49-7561-98343-0  
Fax: +49-7561-98343-699

**SciCan Ltd.**

1440 Don Mills Road  
Toronto, ON M3B 3P9, CANADA  
Telephone: +1-416-445-1600  
Fax: +1-416-445-2727  
Free of charge: +1-800-667-7733  
E-mail: techservice.ca@scican.com

**SciCan Medtech AG**

Alpenstrasse 16  
6300 Zug, SWITZERLAND  
Telephone: +41-41-727-70-27  
Fax: +41-41-727-70-29

**SciCan Inc.**

701 Technology Drive  
Canonsburg, PA 15317, USA  
Telephone: +1-724-820-1600  
Fax: +1-724-820-1479

**MICRO-MEGA® SA**

5-12 rue du Tunnel - BP 1353  
25006 BESANCON Cedex  
FRANCE  
TEL: +33 (0)3 81 54 42 42  
FAX: +33 (0)3 81 54 42 30  
mmb@micro-mega.com

## Table of Contents

1	Important Information.....	5
	Symbols used.....	5
	Safety measures.....	7
	General warnings.....	8
	Warnings for reprocessing (cleaning – disinfection – care – sterilization).....	10
	Intended use – purpose.....	11
	Target group.....	12
2	Technical data.....	13
3	Reprocessing methods.....	18
	Preparation at the location of use.....	18
	Cleaning.....	18
	Sterilization in accordance with ISO 17665.....	20
	Storage.....	23
4	Installation and Commissioning.....	24
	Cooling media supply.....	24
5	Operation.....	25
	Connecting the handpiece.....	26
	Removing handpiece.....	27
	Inserting cutters or diamond grinders.....	28
	Removing cutters or diamond grinders.....	30

## Table of Contents

Inserting rubber cups or polishing brushes (PSI / PSO) .....	31
Removing the rubber cups or polishing brushes (PSI / PSO) .....	32
6 Maintenance and repair .....	33
Care products and maintenance systems.....	33
Care with maintenance systems: <i>STATMATIC</i> / <i>STATMATIC PLUS</i> .....	34
Care with dental care spray .....	35
7 Guarantee.....	36


## 1 Important Information

Dear user,  
we hope that you enjoy using your new quality product. In order to ensure that you use the product correctly, efficiently and safely, please adhere to the following information.

### Symbols used

In the manual/on the product		On the packaging	
	Situation which can result in injury or danger to persons, when the instructions are not followed.		Handle with care
	Important information for users and technicians		Keep dry
	Instructions for disposal		Maximum stacking load
	Steam-sterilizable up to 135 °C (275 °F)		Temperature range
	CE marking		Humidity

## 1 Important Information

In the manual/on the product		On the packaging	
	Thermodisinfestable		Atmospheric pressure
	Manufacturer		Quantity
	Observe instructions for use		
	Follow instructions for use!		

Product subject to technical modifications and design deviations.

No retrofit warranty on existing devices.

The manufacturer does not accept any responsibility for damages caused by:

- external effects (poor quality of media, impact damage, cleaning and care with unsuitable products or faulty installation)
- use of wrong information
- unintended use
- improper repairs

## 1 Important Information

Repair and maintenance, except for the activities described in these instructions for use, must only be performed by qualified professionals.

- Changes made by third parties render the approval null and void.
- Use original parts and spare parts only.

The guarantee becomes ineffective if defects or their consequences are caused by manipulation or modification of the product by the customer or third parties not authorized by us.

Disposal of devices and accessories after use:

Based on EC directive (WEEE 2012/19/EU) on waste electrical and electronic equipment, we hereby inform you that this product is subject to the aforementioned directive and must be disposed of through special channels. Contact the manufacturer or dental supply shop for further information.

### Safety measures

#### **WARNING! Risks when using handpieces with electrical micro-motors!**

Electrical micro-motors generate significantly more energy than conventional air turbines and air motors. Due to the high torques and speeds, poorly serviced, damaged or misused handpieces may lead to overheating which could inflict severe burns on the patient.

▷ Please observe the following.

## 1 Important Information

- ▶ Carefully follow the maintenance instructions.
- ▶ Inspect the handpiece for external damage before each use.
- ▶ Run a test with the handpiece before each use. Especially pay attention to unusual overheating, conspicuous noises and vibrations.
- ▶ If a handpiece is found to be faulty, its use must be discontinued immediately.
- ▶ Never lift cheeks or tongues with the handpiece.

The following people are authorized to perform repair and maintenance works on our products:

- ✓ Our technicians
- ✓ Technicians who have been specially trained by us

### General warnings

- ▶ Ensure that all of the rotation instruments used are in a perfectly-functioning condition.
- ▶ Ensure that the handpiece is securely engaged with the motor coupling and that the rotation instrument is fastened correctly before starting the handpiece.
- ▶ Use a motor that is in a perfectly-functioning condition.
- ▶ Avoid using a handpiece that is unusually warm.

## 1 Important Information

- ▶ It is recommended that appropriate protective gloves and goggles be worn when working with the handpiece and the rotation instruments.
- ▶ Use a rubber dam during patient procedures.
- ▶ During a long treatment period, clean and check the instrument every ten minutes.
- ▶ The handpiece may only be disassembled by an approved retailer or the manufacturer.
- ▶ **Do not** press the push-button of the handpiece during operation.
- ▶ Adjust the rotational speed of the motor to the handpiece used.
- ▶ Do not use rotation instruments (cutter, grinder) which are blunt, bent or damaged.
- ▶ Use SNAO 200 L handpieces exclusively with spray cooling.
- ▶ Should a problem arise with the handpiece, consult our customer service department.
- ▶ Discontinue using the handpiece in the event of damage, overheating, unusual operating noises, excessive vibrations, an excessive lateral run-out of the instrument or an inadequate fixation of the instrument. In order to prevent further damage, consult our customer service department or your local retailer of this product.

## 1 Important Information

### Warnings for reprocessing (cleaning – disinfection – care – sterilization)

- ▶ Clean and sterilize the handpiece prior to the first use and after each patient operation and prior to disposal.
- ▶ Remove the rotation instruments (drills, files, burs) before reprocessing.
- ▶ Only use disinfecting cleaning products that are designed for disinfecting dental handpieces and contra-angles. Do not use setting cleaning agents.
- ▶ Adhere to the concentrations and immersion periods recommended by the manufacturer of the pre-disinfection product.
- ▶ Do not use vinegar to clean the spray heads.
- ▶ Open the sterilization pouches directly before use.
- ▶ Excessive use of the cleaning products (dental care spray, cleaning and disinfection products) may possibly lead to the permeation of these products into the motor and may impair its functionality. The instructions for care for the respective product must be observed.
- ▶ If the handpiece will not be used for a long period of time, clean it as per the instructions after its last use and store it in a dry location.
- ▶ Observe the national guidelines and provisions for the processing of medical products.

## 1 Important Information

### Intended use – purpose

This medical product:

- is only intended for dental treatment in the area of dentistry. Any form of misuse or change of the product is forbidden and may lead to potential dangers.
- is intended for the following applications:
  - SNAO 200 L / 40 ST: The removal of decayed matter, cavity and crown preparations, the removal of fillings and surface finishing of tooth and restoration surfaces.
  - SNAO 40 / 40 L / 10 / 10 L: Cavity preparations, caries excavation, endodontics, surface finishing of tooth and restoration surfaces.
  - SNAO PSI / PSO: Prophylaxis treatment.
- is a medical product intended for use according to the applicable national legal provisions.

## 1 Important Information

According to these regulations this medical product must only be used for the described applications by professional technicians and under consideration of the following issues:

- applicable industrial safety regulations,
- applicable accident prevention measures and
- these instructions for use.

According to these regulations it is the user's responsibility:

- to use only faultless work equipment,
- to ensure the correct purpose of use,
- to protect the patient and others from danger and
- to avoid contamination by the product.

### Target group

This document is directed to the dentist and the dentist's assistant.

The chapter "Commissioning" is additionally directed to the service technician.

## 2 Technical data

	<b>SANAO 200 L</b>	<b>SANAO 40 / 40 L</b>	<b>SANAO 10 / 10 L</b>
Max. drive speed	40,000 min <sup>-1</sup>	40,000 min <sup>-1</sup>	40,000 min <sup>-1</sup>
Transmission	1 : 5	1 : 1	5 : 1
Labeling	Red dots	Blue dots	Green dots
Max. rotation speed	200,000 min <sup>-1</sup>	40,000 min <sup>-1</sup>	8,000 min <sup>-1</sup>
Lighting system	Glass light rod	Glass light rod (only SANAO 40 L)	Glass light rod (only SANAO 10 L)
Chuck	Push-button chuck	Push-button chuck	Push-button chuck
Classification as per 93/42/EEC	Class IIa	Class IIa	Class IIa
Connection of instruments	ISO 3964 / DIN 13940	ISO 3964 / DIN 13940	ISO 3964 / DIN 13940

## 2 Technical data

	<b>SANAO 200 L</b>	<b>SANAO 40 / 40 L</b>	<b>SANAO 10 / 10 L</b>
Tools	Carbide cutter/diamond grinder made of steel or carbide according to ISO 1797-1, type 3	Carbide cutter/diamond drill made of steel or carbide according to ISO 1797-1, type 1	Carbide cutter/diamond drill made of steel or carbide according to ISO 1797-1, type 1
- Shaft diameter	1.59 – 1.60 mm	2.334 – 2.350 mm	2.334 – 2.350 mm
- Max. total length	25 mm	34 mm	34 mm
- Shaft clamping length	min. 11 mm		
- Cutting edge diameter	max. 2 mm		
Air cooling	8 NI/min 2.5 bar (36.3 psi)	8 NI/min 2.5 bar (36.3 psi)	8 NI/min 2.5 bar (36.3 psi)
Spray air	1.8 NI/min 2.7 bar (39.2 psi)	1.8 NI/min 2.7 bar (39.2 psi)	1.8 NI/min 2.7 bar (39.2 psi)
Spray water	50 ml/min 2 bar (29 psi)	50 ml/min 2 bar (29 psi)	50 ml/min 2 bar (29 psi)

## 2 Technical data

	<b>SANAO 40 ST</b>	<b>SANAO PSI / PSO</b>
Max. drive speed	40,000 min <sup>-1</sup>	40,000 min <sup>-1</sup>
Transmission	1 : 1	5 : 1
Labeling	Blue ring	Green dots
Max. rotation speed	40,000 min <sup>-1</sup>	8,000 min <sup>-1</sup>
Classification as per 93/42/EEC	Class IIa	Class IIa
Connection of instruments	ISO 3964 / DIN 13940	ISO 3964 / DIN 13940

## 2 Technical data

	<b>SANAO 40 ST</b>	<b>SANAO PSI / PSO</b>
Tools	<p>Handpiece drill/contra-angle drill made of steel or carbide according to ISO 1797-1, type 2</p> <ul style="list-style-type: none"> <li>- Shaft diameter: 2.334 – 2.350 mm</li> <li>- With drill stop: Shaft clamping length min. 12 mm / max. 22 mm</li> <li>- Without drill stop: Shaft clamping length min. 30 mm / max. 44.5 mm</li> </ul>	<p>Rubber cups/polishing brushes</p> <ul style="list-style-type: none"> <li>- <b>SANAO PSI:</b> With screw-in function as per DIN EN ISO 13295</li> <li>- <b>SANAO PSO:</b> With a snap-on latch as per DIN EN ISO 13295, latch type 5</li> </ul>
Air cooling	8 NI/min      2.5 bar (36.3 psi)	Not applicable
Spray air	1.8 NI/min      2.7 bar (39.2 psi)	Not applicable
Spray water	50 ml/min      2 bar (29 psi)	Not applicable

## 2 Technical data

### Ambient conditions

---

Installation location:	Permitted in interior rooms
Ambient temperature:	10 – 40 °C (50 – 104 °F)
Relative humidity:	30 – 75 %

### Storage and transport conditions

---

Ambient temperature:	Transport: -30 – 70 °C (-22 – 158 °F) Storage: 0 – 40 °C (32 – 104 °F)
Relative humidity:	15 – 93 %
Air pressure:	No restriction
Keep dry!	

Subject to technical changes.

### 3 Reprocessing methods

In the following the reprocessing methods according to ISO 17664 are described.

#### Preparation at the location of use

**WARNING! Danger from non-sterile products!**

Contaminated medical products pose a threat of infection.

- ▷ Take suitable action to protect personnel.
- ▷ Immediately remove residues of cement, composite or blood at the place of use.
- ▷ Transport the medical product in dry condition for reprocessing.
- ▷ Do not soak in solutions.
- ▷ The medical product must be reprocessed right after the treatment.

#### Cleaning

**NOTE! Malfunctions by cleaning in ultrasonic device!**

Defects of the medical product.

- ▷ The medical product must not be placed in a cleaning basin or ultrasonic device.

### 3 Reprocessing methods

#### Option I

#### Manual exterior cleaning

Necessary accessories:

- Potable water 30 °C ± 5 °C (86 °F ± 9 °F)
- Brush, e.g., medium-hard toothbrush

Brush under running potable water for 20 seconds.

#### Manual interior cleaning and maintenance (for U.S. market)

Use *STATCARE* spray in accordance with the instructions for use.

#### Option II

#### Automated internal and external cleaning and disinfection in a washer-disinfector according to ISO 15883

#### NOTE! Wet medical product after reprocessing!

Defects with the use of wet medical product.

- ▷ Choose suitable cleaning and drying method.
- ▷ Ensure that the medical product is completely dry after reprocessing.

In some countries the applicable national guidelines require automated reprocessing in a validated washer-disinfector (WD) such as, for instance, SciCan HYDRIM M2 or C61! Observe instructions for use.

#### Validated parameters for thermal disinfection

Detergent:	TPH 5949 (S&M)
PH-Value:	max. 10
Program:	Thermal disinfection
Temperature:	≥ 93 °C
Exposure time:	≥ 5 minutes
Drying time:	≥ 15 minutes

### 3 Reprocessing methods

#### Option I

##### **Automated interior cleaning and maintenance**

Use an automated cleaning and maintenance device, e.g., SciCan *STATMATIC PLUS*, in accordance with the instructions for use.

##### **Disinfection**

The medical product can be disinfected in a steam sterilizer, e.g., SciCan *STAT/M*.

- ▷ Observe national and regional guidelines.

##### **Sterilization in accordance with ISO 17665**

##### **Packaging**

The packaging of sterile goods must be adapted to the instrument in order to have a tension-free packaging.

The packaging of sterile goods must fulfil the applicable standards regarding quality and use and must be suitable for the sterilization method!

- ▶ Package the medical products individually in packaging for sterile goods.
- ▶ Perform sterilization in a steam sterilizer (autoclave) in accordance with EN 13060 following ISO 17665-1.

#### Option II

- ▷ Use only WDs with integrated drier.
- ▷ Select an appropriate program on the WD with drying cycle.
- ▷ Observe national guidelines.

##### **Automated maintenance**

Use an automated maintenance device, e.g., SciCan *STATMATIC*, in accordance with the instructions for use.

### 3 Reprocessing methods

**NOTE! Premature wear and malfunctions by improper maintenance and care!**

Reduced life span of the product.

▷ Treat medical product with dental cleaning products before each sterilization cycle.

**NOTE! Contact corrosion caused by moisture!**

Damage to the medical product.

▷ Immediately remove the medical product from the steam sterilizer after the sterilization cycle.

The medical product can withstand a temperature of up to 138 °C (280.4 °F).

The medical product can be sterilized in steam sterilizers.

Observe national and regional guidelines.

Observe the instructions of the autoclave manufacturer.

Sterilization in a standards-compliant (e.g. EN 13060) steam sterilizer with validated B or S cycle.

### 3 Reprocessing methods

Recommended sterilization methods or parameters are as follows:

Pressure impulse method S cycle	134 °C (273.2 °F) Minimum pressure 300 kPa	Minimum sterilization time	3.5 - 5 - 18 minutes*	Minimum drying time	15 minutes
Vacuum method B cycle	134 °C (273.2 °F) Minimum pressure 300 kPa	Minimum sterilization time	3.5 - 5 - 18 minutes*	Minimum drying time	10 minutes

\* Depending on the type of use. National guidelines have to be observed.

Follow the instructions of the autoclave manufacturer regarding the sterilization of this medical products for autoclaves with S cycle method.

Pursuant to U.S. and Canadian standards:

Dynamic Air Removal: Unwrapped	132 °C (270 °F)	Minimum sterilization time = 3 minutes	Drying time not applicable
	134 °C (273 °F)	Minimum sterilization time = 3.5 minutes	Drying time not applicable
	135 °C (275 °F)	Minimum sterilization time = 3 minutes	Drying time not applicable
Dynamic Air Removal: Wrapped or pouched	132 °C (270 °F)	Minimum sterilization time = 4 minutes	Drying time = 60 minutes
	134 °C (273 °F)	Minimum sterilization time = 3.5 minutes	Drying time = 60 minutes
	135 °C (275 °F)	Minimum sterilization time = 3 minutes	Drying time = 60 minutes

**NOTE! This table represents the variation in sterilizer manufacturers' recommendations for exposure at different temperatures!**

▷ For a specific sterilizer, consult only that manufacturer's recommendations.

### 3 Reprocessing methods

#### Storage

Store reprocessed products in a dust-protected, dry, dark and cool room which is as sterile as possible.

Please observe the use-by date of the sterile product!

## 4 Installation and Commissioning

### **WARNING! Risks posed by unsterilized products!**

Risk of infection to user and patient.

- ▷ Before putting the device into operation for the first time and after each use, reprocess the medical product accordingly and sterilize if necessary.

### **NOTE! Malfunction of the handpiece as a result of a low ambient temperature!**

Possible malfunctions of the handpiece.

- ▷ Before putting into operation, return cooled products to a temperature of 20 – 25 °C (68 – 77 °F).

### **NOTE! Damage may occur as a result of polluted and moist cooling air!**

Polluted and moist cooling air may lead to malfunctions and premature bearing wear.

- ▷ Ensure that the cooling air is dry, clean and unpolluted as per ISO 7494-2 standards.

## Cooling media supply

### **CAUTION! Overheating of tooth as a result of a low volume of water!**

Thermal damage of the pulp.

- ▷ Set water volume to at least 50 ml/min for the spray cooling.
- ▷ Inspect spray water channels and clean spray nozzles if necessary.

## 5 Operation

**WARNING! Dangers will result from damage, unusual operating noises, excessively strong vibrations, excessive heating or if the cutter or grinder are not tightly fastened!**

Dangers to the user and patient.

▷ Discontinue use and inform the service department.

**WARNING! Risk of burns from the hot head or lid of the instrument!**

If the instrument overheats, burns may occur in the mouth area.

▷ Never touch soft tissue using the head or lid of the instrument.

**CAUTION! Risk posed by an incorrectly stored instrument!**

Clamped cutters or grinders may lead to injury or infection.

Damage to the chuck may occur if the instrument is dropped.

▷ Correctly return instrument to the tray without the cutter or grinder after the treatment.

**NOTE! Premature wear and malfunctions may arise as a result of incorrect storage before long periods of non-use!**

Shorter product life cycle.

▷ Before long periods of non-use, the medical product should be cleaned according to the instructions, treated and stored in a dry location.

## 5 Operation

### Connecting the handpiece

**NOTE! The handpiece may become loose during treatment!**

An incorrectly-latched handpiece may detach from the motor coupling and fall down.

▷ Before each treatment, check whether the handpiece is securely latched with the motor coupling by gently pulling it.

**NOTE! Removing and returning the handpiece while the drive motor is rotating.**

Damage to the catch.

▷ Never remove or return the handpiece while the drive motor is rotating.

**NOTE! When used together with a drive motor (SANA O 40 ST)!**

The handpiece may become blocked.

▷ Only use handpiece when the chuck is in a closed position.


**NOTE! Damage (SANA O 200 L / 40 / 40 L / 10 / 10 L)!**

Pressing the foot switch while connecting and removing the handpiece may lead to damage to the handpiece and to the motor coupling.


▷ Do not connect or remove the handpiece while pressing the foot switch.

## 5 Operation

1. Slightly moisten o-rings on the motor coupling using a dental care spray.
2. Insert the handpiece in the motor coupling and rotate until the latch audibly engages.
3. Ensure that the handpiece is securely inserted in the coupling by pulling.

200 L / 40 / 40 L / 10 / 10 L / PSI / PSO

40 ST

### Removing handpiece

Unlatch handpiece from motor coupling by rotating slightly and remove in the direction of the axis.

## 5 Operation

### Inserting cutters or diamond grinders

Only use carbide cutters or diamond grinders that meet the criteria specified under chapter 2 "Technical Data"!

#### **WARNING! Use of unauthorized cutters or grinders!**

This may cause injury to the patient or damage to the handpiece.

- ▷ Observe the instructions for use and only use cutters or grinders according to the intended use.
- ▷ Only use cutters or grinders which do not deviate from the specified data.

#### **WARNING! Risk of injury from cutters or grinders!**

Infections or cutting injuries.

- ▷ Wear protective gloves or finger guards.

#### **CAUTION! Injuries when using worn cutters or grinders!**

Cutters or grinders could fall out and injure the patient during the treatment.

- ▷ Never use cutters or grinders with worn shafts.

## 5 Operation


### **CAUTION! Risks posed by faulty chuck!**

The cutter or grinder could fall out and lead to injuries.

- ▷ When inspecting, inserting and removing, wear protective gloves or finger guards in order to prevent risk of injury or infection.
- ▷ Pull on the cutter or grinder to check whether the chuck is functioning correctly and to ensure that the cutter or grinder is in a fixed position.


How to use the handpieces: **200 L / 40 / 40 L / 10 / 10 L**

1. Firmly press the push-button using your thumb and simultaneously insert the cutter or grinder as far as it will go.
2. Ensure that the cutter or grinder is in a fixed position by pulling it.
3. Ensure that the push-button lies evenly and is not tilted.

How to use the handpieces: **40 ST**

1. Turn the clamping ring anti-clockwise and insert the cutter or grinder into the chuck as far as it will go.
2. Re-tighten the clamping ring by turning it clockwise.
3. Ensure that the cutter or grinder is in a fixed position by pulling it.

## 5 Operation

### Removing cutters or diamond grinders

#### **CAUTION! Risk of injury from the rotating cutter or grinder!**

Cut injuries and damage to the chuck.

- ▷ Do not touch rotating cutters or grinders.
- ▷ Never press the push-button while the cutter or grinder is rotating.
- ▷ Remove cutter/grinder from the handpiece after the treatment has ended in order to prevent injury and infections when returning to the holder.


How to use the handpieces: **200 L / 40 / 40 L / 10 / 10 L**

- ▶ After the cutter or grinder has been switched off, firmly press the push-button using your thumb and simultaneously pull the cutter or grinder out.

How to use the handpieces: **40 ST**

1. After the cutter or grinder has been switched off, turn the clamping ring anti-clockwise and pull the cutter or grinder out.
2. Re-tighten the clamping ring by turning it clockwise.

## 5 Operation

### Inserting rubber cups or polishing brushes (PSI / PSO)

! Only use rubber cups or polishing brushes which meet those criteria specified under chapter 2 "Technical Data"!

**WARNING! Use of unauthorized inserts!**

This may cause injury to the patient or damage to the handpiece.

▷ Observe the instructions for use and the intended purpose of the insert.

**WARNING! Risk of injury from the insert!**

Infections or cutting injuries.

▷ Wear protective gloves or finger guards.

1. Block the catch.
2. Screw the insert into the head (**SANAO PSI**) or slide it into place (**SANAO PSO**).
3. Check that the rubber cup and polishing brush are in a fixed position by pulling.

## 5 Operation

### Removing the rubber cups or polishing brushes (PSI / PSO)

**WARNING! Risk posed by the rotating insert!**

Cutting injuries.

- ▷ Do not touch the rotating insert.
- ▷ Remove insert from the handpiece after the treatment has ended in order to prevent injury and infections when returning to the holder.

- ▶ After switching off, unscrew the insert (SANA O PSI) or take it off (SANA O PSO).

## 6 Maintenance and repair

### **CAUTION! Risk of injuries and damage to material by improper maintenance!**

- ▷ Have maintenance work only performed by an authorized dealer.
- ▷ Have the device sent to maintenance once a year after the guarantee period has elapsed or after 250 processing cycles.

To guarantee a flawless function of the medical product:

- ▶ Treat the medical product in accordance with the processing methods described in the instructions of use and use the care products and maintenance system described therein.
- ▶ Determine an internal maintenance interval for the practice, at which the medical product is evaluated by an authorized dealer in terms of cleaning, maintenance and function. This maintenance interval depends on the frequency of use and must be adapted accordingly.

### **Care products and maintenance systems**

### **CAUTION! Sharp milling cutter or grinder in medical product!**

Risk of injury by sharp and/or pointy milling cutter or grinder.

- ▷ Remove milling cutter or grinder.

### **NOTE! Premature wear and malfunctions by improper maintenance and care!**

Reduced life span of the product.

- ▷ Regularly perform proper maintenance.

## 6 Maintenance and repair

### Care with maintenance systems: *STATMATIC* / *STATMATIC PLUS*

The medical product must be maintained after each cleaning and before each sterilization.

1. Remove milling cutter or grinder.
2. Maintain medical product in accordance with instructions for use.

#### Collet care

Perform collet care once a week, see instructions for use "*STATMATIC* / *STATMATIC PLUS*", Chapter "Collet care program".

## 6 Maintenance and repair

### Care with dental care spray

The medical product must be maintained after each cleaning and before each sterilization.

1. Remove milling cutter or grinder.
2. Cover medical product with maintenance bag.
3. Place medical product on the cannula and activate spray button for one second.

### Collet care

Perform collet care once a week.

- ▶ Remove milling cutter or grinder and spray with the tip of the spray nozzle into the opening.

## 7 Guarantee

We offer a guarantee to the user covering all functional, material and manufacturing defects. This guarantee is valid for a period of 2 years from invoice date (the guarantee period may vary in your country, please see your purchasing documents). In the case of justifiable claims we or an authorized representative shall fulfil the guarantee obligations of the company by repairing or replacing the product free of charge. All other claims of any kind, in particular compensation for damages and interest receivable, are excluded from this guarantee.

We hold no responsibility for damage to materials or persons or their consequences, which are caused by the following:

- excessive wear,
- improper handling,
- non-compliance with the instructions of installation, use and maintenance,
- unusual chemical, electrical or electrolytic influences and
- improper connection to air, water and current supply. The guarantee does not extend to flexible glass fibre connectors or parts made of plastic.

This guarantee can become null and void, if damages and their consequences are caused by the improper use of the product or modifications, which were made by persons unauthorized by us. Claims from this guarantee are only taken in consideration, when invoice or delivery slip are submitted together with the product, which clearly indicate the purchase date as well as the product designation and serial number.

(DE = original)

# SANAO

## Handstück

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EN  
DE  
FR

CS / IT / NL / ...

via: [www.scican.uk.com/downloads](http://www.scican.uk.com/downloads)

200 L  
40 / 40 L  
10 / 10 L  
PSO  
PSI  
40 ST

Rx only

# SciCan

Your Infection Control Specialist™

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