

## Safety Data Sheet

Revision Date: 08.05.2015

Version: 0

Date of Print: 28.05.2015

<b>1. Product identification/preparation/company undertaking</b>
<b>Product details</b> <b>Trade name: PolySil Light ( Fast and Reg. Set)</b> <b>Use or preparation:</b> Impression material <b>Supplier:</b> SciCan Ltd 1440 Don Mills Road Toronto, Ontario, Canada M3B3P9 <b>Emergency information:</b> 800-667-7733
<b>2. Hazards identification</b>
<b>Classification of the substance or mixture</b> <b>Classification according to Regulation (EC) 1272/2008:</b> void <b>Classification according to Directive 67/548/EEC or Directive 1999/45/EC:</b> void <b>Information concerning particular hazards for human and environment:</b> The product does not have to be labeled due to the calculation procedure of the „General classification guideline for preparations of the EU“ in the latest valid version. <b>Label elements:</b> <b>Labeling according to EU guidelines:</b> Medical devices as defined in Directive 93/42/EEC and which are invasive or used in direct physical contact with the human body; are exempted from the provisions of dangerous substance Directive 1999/45/EWG. <b>Other hazards</b> <b>Results of PBT and vPvB assessment</b> <b>PBT:</b> not applicable <b>vPvB:</b> not applicable
<b>3. Composition/Information of Ingredients</b>
<b>Chemical characterization:</b> <b>Description:</b> mixture of polydimethyl polymethyl vinyl siloxane, polydimethyl polymethyl hydrogen siloxane and silica <b>Hazardous components:</b> none
<b>4. First-aid measures:</b>
<b>General information:</b> no special measures required <b>after skin contact:</b> wash with plenty of water and soap <b>after eye contact:</b> rinse with plenty of water and consult a doctor <b>after swallowing:</b> If symptoms persist consult a doctor
<b>5. Fire fighting measures</b>
<b>Extinguishing media:</b> CO <sub>2</sub> , water <b>Protective equipment:</b> no special measures required

## Safety Data Sheet

Revision Date: 08.05.2015  
Date of Print: 28.05.2015

Version: 0

<b>6. Accidental release measures</b>			
<b>Personal precautions:</b> not required			
<b>Environmental precautions:</b> no special measures required			
<b>Methods for cleaning up:</b> pick up mechanically			
<b>Additional information:</b> none			
<b>7. Handling and storage</b>			
<b>Handling:</b>			
<b>Information for safe handling:</b> For dental use only. No special measures required.			
<b>Recommendation for fire and explosion protection:</b> No special measures required			
<b>Storage:</b>			
<b>Requirements at storerooms and containers:</b> No special measures required			
<b>Requirements for storage with other products:</b> not required			
<b>Further information on storage conditions:</b> none			
<b>8. Exposure controls/personal protection</b>			
<b>Components with limits of values to be supervised at the workplace:</b>			
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace			
<b>Additional information:</b> The lists valid during the making were used as basis.			
<b>Personal protective equipment:</b>			
<b>General measures of protection and hygiene:</b> normal hygienic measures			
<b>Respiratory protection:</b> not required			
<b>Protection of hands:</b> not required			
<b>Eye protection:</b> not absolute required			
<b>9. Physical and chemical properties</b>			
<b>Appearance:</b>			
<b>Form:</b>	paste		
<b>Color:</b>	different		
<b>Odor:</b>	odorless		
	<u>Value</u>	<u>Unit</u>	<u>Method</u>
<b>Information on change in the physical state</b>			
<b>Melting point/melting range:</b>	n.a.		°C
<b>Boiling point/boiling range:</b>	n.a.		°C
<b>Flash point:</b>	n.a.		°C
<b>Autoignition temperature:</b>	n.a.		
<b>Danger of explosion:</b>	none		
<b>Density:</b>	1,24 (20°C)	g/cm <sup>3</sup>	
<b>Vapor pressure:</b>	n.a.	mbar	
<b>Viscosity:</b>	paste		
<b>pH:</b>	neutral		
<b>Solubility in/miscibility with:</b>	Partially soluble in toluene, petrol ether		
<b>Water:</b>	insoluble		
<b>Content of solvents:</b>	none		
<b>Organic solvents:</b>	----		
<b>Water:</b>	----		
<b>Content of solids:</b>	n.a.		

## Safety Data Sheet

Revision Date: 08.05.2015

Version: 0

Date of Print: 28.05.2015

<p><b>10. Stability and reactivity</b></p> <p><b>Conditions to avoid:</b> No decomposition if used according to specification</p> <p><b>Hazardous decomposition products:</b> none under normal conditions of storage and use</p>
<p><b>11. Toxicological information</b></p> <p><b>Acute toxicity:</b></p> <p><b>Primary irritation:</b> <b>skin:</b> no irritating effect <b>eye:</b> no irritating effect</p> <p><b>Additional toxicological information:</b> When used and handled according to specifications, the product does not have any harmful effects to our experience and the information provided to us.</p>
<p><b>12. Ecological information</b></p> <p><b>General information:</b> Avoid transfer into environment</p> <p><b>Classification of water endangerment:</b> WGK 1 (German regulation) slightly hazardous for water</p>
<p><b>13. Disposal consideration</b></p> <p><b>Product:</b></p> <p><b>Recommendation:</b> Small quantities can be disposed of with household waste. Use proper landfill disposal or incineration in accordance with local, state and federal regulations.</p> <p><b>Uncleaned packaging:</b></p> <p><b>Recommendation:</b> Disposal must be made according to official regulations.</p>
<p><b>14. Transport information</b></p> <p><b>Land transport ADR, RID: -</b> <b>Maritime transport: IMDG-Code: -</b> <b>Air transport ICAO-TI/IATA-DGR: -</b></p>
<p><b>15. Regulatory information</b></p> <p><b>Classification according to EC-guidelines:</b> The product is a medical device according to EC-directive 93/42 EEC.</p> <p><b>Chemical safety assessment:</b> A Chemical Safety Assessment has not been carried out.</p>
<p><b>16. Other information:</b></p> <p>The contents of this document have been prepared in accordance with Regulation (EC) No 1907/2006 (REACH), Annex II (453/2010). The above information is based on our present day knowledge and relates solely to the safety requirements of the product. The data do not signify any warranty with regards to products properties. However users of the product should satisfy themselves that the information given is sufficient and correct for their specific circumstances of use.</p> <p>Reference: PolySil Light ( Fast and Reg. Set)- EU &amp; US SDS; Date of issue: 30.10.2013; Preparation by Dr. Jürgen Engelbrecht, Tel: 04121/483-0</p>