

AQUASTAT™ – Perfect for your office needs.

Every drop. Perfect, steam distilled water.



- **Easy-to-use and energy efficient**

Simply plug it in. No special wiring or plumbing required. Process 3.78 litres (1 gallon) of distilled water in only 4.25 hours. Lightweight design includes conveniently removeable boiling chamber for easy maintenance.

- **Safe, clean water from any source**

Whether your water source is spring fed, municipally piped, soft or hard, the AQUASTAT assures you of fresh, safe, clean, steam-distilled water every time.

- **Protect your equipment and remove water contaminants**

Vaporizes tap water to kill and remove 99% of most contaminants that could damage your equipment.



AQUASTAT

AQUASTAT steam distiller produces perfect, steam distilled water for your autoclave – ideal for STATIM®

Save the environment, stop wasting money on expensive bottled water and avoid damaging equipment from poor quality distilled water.

Weight: 4.5 kg/10 lbs
Dimensions (L x W x H):
 406 mm x 245 mm x 384 mm
 16" x 9 5/8" x 15" 1/8"

Power: 230 VAC, 50Hz, 750W
P/N: W10241S (EU Plug)
 W10243S (Swiss Plug)
 W10244S (UK Plug)

AQUASTAT™ – Ideal für Ihre Praxis.

Perfekt dampf-destilliertes Wasser – jeder Tropfen.



- Bedienerfreundlich und energieeffizient**
 Einfach einstecken. Keine spezielle Verkabelung oder Sanitärinstallation erforderlich. Produziert 3,78 Liter destilliertes Wasser in nur 4,25 Stunden. Das leichtgewichtige Gerät verfügt - zur einfacheren Handhabung - über einen entnehmbaren Siedebehälter.
- Sauberes Wasser ganz gleich woher**
 Egal ob Sie Quellwasser, städtisches Leitungswasser, weiches oder hartes Wasser verwenden, mit dem AQUASTAT haben Sie immer frisches und sauberes dampf-destilliertes Wasser.
- Schützen Sie Ihre Geräte und entfernen Sie Schadstoffe im Wasser**
 Verdampft Leitungswasser und entfernt dabei 99% der häufigsten Verunreinigungen die Ihre Geräte beschädigen können.



AQUASTAT

Das AQUASTAT Destilliergerät produziert perfekt dampf-destilliertes Wasser für Ihren Autoklav – ideal für STATIM®

Schützen Sie die Umwelt, geben Sie kein Geld mehr für teures abgefülltes Wasser aus und schützen Sie Ihre Geräte vor Schäden durch destilliertes Wasser von schlechter Qualität.

Gewicht: 4,5 kg
Maße (L x B x H):
 406 mm x 245 mm x 384 mm

Anschlussdaten: 230 V, 50Hz, 750W
Art.-Nr.: W10241S (EU)
 W10243S (CH)
 W10244S (UK)